


## ORDER SAMPLES

Visit our website to request a sample of our surgical instruments and experience firsthand the commitment to quality that defines Makvin Surgical.




# MakEdge

Single Use  
Ophthalmic  
Micro Surgical Knives



## CONTACT US



## Makvin Surgical LLC

Phone : +1 423 750 4247  
E-mail : [sales@makvinsurgical.com](mailto:sales@makvinsurgical.com)  
Website : [www.makvinsurgical.com](http://www.makvinsurgical.com)

**MAKVIN**  
Surgical



## Who We Are

---

Makvin Surgical is an US based company, supplies single use surgical instruments designed for use by healthcare professional. We emphasize reliability and consistency across our products.

### Mission

---

We strive to support surgeons,  
in their pursuit of better quality of care.

### Vision

---

We aim to be a global leader in surgical products,  
trusted for consistency and quality.

### Values

---

#### Integrity

We do what's right

#### Accountability

We do what we say

#### Empathy

We listen, understand and respond

#### Transparency

We communicate clear and honest

## Why MakEdge?

---

MakEdge ophthalmic surgical knives are single-use instruments intended for use by surgeons and trained healthcare professionals. Manufactured using automated electro-etching and electro-polishing processes to support uniform blade formation and a refined edge finish, each knife features a sharp, clean cutting edge and a matte finish that helps reduce glare under operating room lighting—supporting a clearer view of the surgical field. An ergonomic handle with a textured surface supports a secure grip for confident handling during procedures.

### At a Glance

---

- Single-use
- Sharp, clean cutting edge
- Matte anti-glare finish
- Ergonomic, textured grip
- Depth indicator (selected models)
- Sterile (EO)



# Slit Knives



Single Bevel:

Qty : 6 Pieces/Box

SKU	DESCRIPTION
1418120	1.2MM, ANGLED, SINGLE BEVEL, SLIT KNIFE
1418180	1.8MM, ANGLED, SINGLE BEVEL, SLIT KNIFE
1418200	2.0MM, ANGLED, SINGLE BEVEL, SLIT KNIFE
1418220	2.2MM, ANGLED, SINGLE BEVEL SLIT KNIFE
1418240	2.4MM, ANGLED, SINGLE BEVEL, SLIT KNIFE
1418250	2.5MM, ANGLED, SINGLE BEVEL, SLIT KNIFE
1418260	2.6MM, ANGLED, SINGLE BEVEL, SLIT KNIFE
1418265	2.65MM, ANGLED, SINGLE BEVEL, SLIT KNIFE
1418275	2.75MM, ANGLED, SINGLE BEVEL, SLIT KNIFE
1418280	2.8MM, ANGLED, SINGLE BEVEL, SLIT KNIFE
1418300	3.0MM, ANGLED, SINGLE BEVEL, SLIT KNIFE
1418320	3.2MM, ANGLED, SINGLE BEVEL, SLIT KNIFE

Double Bevel:

Qty : 6 Pieces/Box

SKU	DESCRIPTION
1428200	2.0MM, ANGLED, DOUBLE BEVEL, SLIT KNIFE
1428220	2.2MM, ANGLED, DOUBLE BEVEL, SLIT KNIFE
1428240	2.4MM, ANGLED, DOUBLE BEVEL, SLIT KNIFE
1428250	2.5MM, ANGLED, DOUBLE BEVEL, SLIT KNIFE
1428260	2.6MM, ANGLED, DOUBLE BEVEL, SLIT KNIFE
1428275	2.75MM, ANGLED, DOUBLE BEVEL, SLIT KNIFE
1428280	2.8MM, ANGLED, DOUBLE BEVEL, SLIT KNIFE
1428285	2.85MM, ANGLED, DOUBLE BEVEL, SLIT KNIFE
1428300	3.0MM, ANGLED, DOUBLE BEVEL, SLIT KNIFE
1428320	3.2MM, ANGLED, DOUBLE BEVEL, SLIT KNIFE

# Clear Cornea Knife



Single Bevel:

Qty : 6 Pieces/Box

SKU	DESCRIPTION
1518220	CLEAR CORNEA KNIFE 2.2MM ANGLED SINGLE BEVEL
1518240	CLEAR CORNEA KNIFE 2.4MM ANGLED SINGLE BEVEL
1518250	CLEAR CORNEA KNIFE 2.5MM ANGLED SINGLE BEVEL
1518265	CLEAR CORNEA KNIFE 2.65MM ANGLED SINGLE BEVEL
1518275	CLEAR CORNEA KNIFE 2.75MM ANGLED SINGLE BEVEL
1518280	CLEAR CORNEA KNIFE 2.8MM ANGLED SINGLE BEVEL
1518300	CLEAR CORNEA KNIFE 3.0MM ANGLED SINGLE BEVEL

Double Bevel:

Qty : 6 Pieces/Box

SKU	DESCRIPTION
1528220	CLEAR CORNEA KNIFE 2.2MM ANGLED DOUBLE BEVEL
1528240	CLEAR CORNEA KNIFE 2.4MM ANGLED DOUBLE BEVEL
1528250	CLEAR CORNEA KNIFE 2.5MM ANGLED DOUBLE BEVEL
1528265	CLEAR CORNEA KNIFE 2.65MM ANGLED DOUBLE BEVEL
1528275	CLEAR CORNEA KNIFE 2.75MM ANGLED DOUBLE BEVEL
1528280	CLEAR CORNEA KNIFE 2.8MM ANGLED DOUBLE BEVEL
1528300	CLEAR CORNEA KNIFE 3.0MM ANGLED DOUBLE BEVEL

## Crescent Knives



Qty : 6 Pieces/Box

SKU	DESCRIPTION
1617200	CRESENT KNIFE 2.0MM STRAIGHT
1617260	CRESENT KNIFE 2.6MM STRAIGHT
1618200	CRESENT KNIFE 2.0MM BEVEL UP
1618220	CRESENT KNIFE 2.2MM BEVEL UP
1618250	CRESENT KNIFE 2.5MM BEVEL UP
1618260	CRESENT KNIFE 2.6MM BEVEL UP

## MVR Knives



Qty : 6 Pieces/Box

SKU	DESCRIPTION
1828190	MVR KNIFE 19 GAUGE ANGLED
1827190	MVR KNIFE 19 GAUGE LANCE
1828200	MVR KNIFE 20 GAUGE ANGLED
1827200	MVR KNIFE 20 GAUGE LANCE
1828210	MVR KNIFE 21 GAUGE ANGLED
1827210	MVR KNIFE 21 GAUGE LANCE
1828230	MVR KNIFE 23 GAUGE ANGLED
1827230	MVR KNIFE 23 GAUGE LANCE
1828250	MVR KNIFE 25 GAUGE ANGLED
1827250	MVR KNIFE 25 GAUGE LANCE


## Stab Knives



Qty : 6 Pieces/Box

SKU	DESCRIPTION
1727150	STAB KNIFE 15 DEGREE
1724150	STAB KNIFE 15 DEGREE ANGLED LEFT
1725150	STAB KNIFE 15 DEGREE ANGLED RIGHT
1727225	STAB KNIFE 22.5 DEGREE
1727300	STAB KNIFE 30 DEGREE
1727450	STAB KNIFE 45 DEGREE

## Handling Instructions



**Taking out the knife:** Take special care to prevent the blade from coming into contact with any part of the packaging. Gently pick up the device and lift the knife from the locking position.

**Replacing the knife:** Hold the knife with the blade edge facing downward. Align it properly with the locking position in the tray. Press down gently until the handle clicks into place, ensuring the knife is securely locked.