

# MakEdge



**Single Use  
Ophthalmic  
Micro Surgical Knives**

# Who We Are

---

Makvin Surgical is a US based company, supplies single use surgical instruments designed for use by healthcare professionals. We emphasize reliability and consistency across our products.

## Mission

---

We strive to support surgeons,  
in their pursuit of better quality of care.

## Vision

---

We aim to be a global leader in surgical products,  
trusted for consistency and quality.

## Values

---

### **Integrity**

We do what's right

### **Accountability**

We do what we say

### **Empathy**

We listen, understand and respond

### **Transparency**

We communicate clear and honest

# Why MakEdge?

MakEdge ophthalmic surgical knives are single-use instruments intended for use by surgeons and trained healthcare professionals. Manufactured using automated electro-etching and electro-polishing processes to support uniform blade formation and a refined edge finish, each knife features a sharp, even cutting edge and a matte finish that helps reduce glare under operating room lighting—supporting a clearer view of the surgical field. An ergonomic handle with a textured surface supports a secure grip for confident handling during procedures.

## At a Glance

---

- Single-use
- Sharp, even cutting edge
- Matte anti-glare finish
- Ergonomic, textured grip
- Depth indicator (selected models)
- Sterile (EO)



## Slit Knives



Depth Indicator



Angle



### Single Bevel:

Qty : 6 Pieces/Box

SKU	DESCRIPTION
KR 12	SLIT KNIFE 1.2MM ANGLED SINGLE BEVEL
KR 18	SLIT KNIFE 1.8MM ANGLED SINGLE BEVEL
KR 20	SLIT KNIFE 2.0MM ANGLED SINGLE BEVEL
KR 22	SLIT KNIFE 2.2MM ANGLED SINGLE BEVEL
KR 24	SLIT KNIFE 2.4MM ANGLED SINGLE BEVEL
KR 25	SLIT KNIFE 2.5MM ANGLED SINGLE BEVEL
KR 26	SLIT KNIFE 2.6MM ANGLED SINGLE BEVEL
KR 265	SLIT KNIFE 2.65MM ANGLED SINGLE BEVEL
KR 275	SLIT KNIFE 2.75MM ANGLED SINGLE BEVEL
KR 28	SLIT KNIFE 2.8MM ANGLED SINGLE BEVEL
KR 30	SLIT KNIFE 3.0MM ANGLED SINGLE BEVEL
KR 32	SLIT KNIFE 3.2MM ANGLED SINGLE BEVEL

(Note : All models of slit knife come with depth indicator)

### Double Bevel:

Qty : 6 Pieces/Box

SKU	DESCRIPTION
KR DB 20	SLIT KNIFE 2.0MM ANGLED DOUBLE BEVEL
KR DB 22	SLIT KNIFE 2.2MM ANGLED DOUBLE BEVEL
KR DB 24	SLIT KNIFE 2.4MM ANGLED DOUBLE BEVEL
KR DB 25	SLIT KNIFE 2.5MM ANGLED DOUBLE BEVEL
KR DB 26	SLIT KNIFE 2.6MM ANGLED DOUBLE BEVEL
KR DB 275	SLIT KNIFE 2.75MM ANGLED DOUBLE BEVEL
KR DB 28	SLIT KNIFE 2.8MM ANGLED DOUBLE BEVEL
KR DB 285	SLIT KNIFE 2.85MM ANGLED DOUBLE BEVEL
KR DB 30	SLIT KNIFE 3.0MM ANGLED DOUBLE BEVEL
KR DB 32	SLIT KNIFE 3.2MM ANGLED DOUBLE BEVEL

(Note : All models of slit knife come with depth indicator)

## Clear Cornea Knives



Depth Indicator



Angle



### Single Bevel:

Qty : 6 Pieces/Box

SKU	DESCRIPTION
CC 22	CLEAR CORNEA KNIFE 2.2MM ANGLED SINGLE BEVEL
CC 24	CLEAR CORNEA KNIFE 2.4MM ANGLED SINGLE BEVEL
CC 25	CLEAR CORNEA KNIFE 2.5MM ANGLED SINGLE BEVEL
CC 265	CLEAR CORNEA KNIFE 2.65MM ANGLED SINGLE BEVEL
CC 275	CLEAR CORNEA KNIFE 2.75MM ANGLED SINGLE BEVEL
CC 28	CLEAR CORNEA KNIFE 2.8MM ANGLED SINGLE BEVEL
CC 30	CLEAR CORNEA KNIFE 3.0MM ANGLED SINGLE BEVEL

(Note : All models of clear cornea knife come with depth indicator)


### Double Bevel:

Qty : 6 Pieces/Box

SKU	DESCRIPTION
CC DB 22	CLEAR CORNEA KNIFE 2.2MM ANGLED DOUBLE BEVEL
CC DB 24	CLEAR CORNEA KNIFE 2.4MM ANGLED DOUBLE BEVEL
CC DB 25	CLEAR CORNEA KNIFE 2.5MM ANGLED DOUBLE BEVEL
CC DB 265	CLEAR CORNEA KNIFE 2.65MM ANGLED DOUBLE BEVEL
CC DB 275	CLEAR CORNEA KNIFE 2.75MM ANGLED DOUBLE BEVEL
CC DB 28	CLEAR CORNEA KNIFE 2.8MM ANGLED DOUBLE BEVEL
CC DB 30	CLEAR CORNEA KNIFE 3.0MM ANGLED DOUBLE BEVEL

(Note : All models of clear cornea knife come with depth indicator)

## Crescent Knives



Qty : 6 Pieces/Box

SKU	DESCRIPTION
CRS 20	CRESCENT KNIFE 2.0MM STRAIGHT BEVEL UP
CRS 26	CRESCENT KNIFE 2.6MM STRAIGHT BEVEL UP
CRU 20	CRESCENT KNIFE 2.0MM ANGLED BEVEL UP
CRU 22	CRESCENT KNIFE 2.2MM ANGLED BEVEL UP
CRU 25	CRESCENT KNIFE 2.5MM ANGLED BEVEL UP
CRU 26	CRESCENT KNIFE 2.6MM ANGLED BEVEL UP


## Stab Knives



Qty : 6 Pieces/Box

SKU	DESCRIPTION
SP 15.0	STAB KNIFE 15 DEGREE
SP 15.0-A-L	STAB KNIFE 15 DEGREE ANGLED (45°) LEFT
SP 15.0-A-R	STAB KNIFE 15 DEGREE ANGLED (45°) RIGHT
SP 22.5	STAB KNIFE 22.5 DEGREE
SP 30.0	STAB KNIFE 30 DEGREE
SP 45.0	STAB KNIFE 45 DEGREE


# MVR Knives



Qty : 6 Pieces/Box

SKU	DESCRIPTION
MVR 19A	MVR KNIFE 19 GAUGE ANGLED
MVR 19	MVR KNIFE 19 GAUGE LANCE
MVR 20A	MVR KNIFE 20 GAUGE ANGLED
MVR 20	MVR KNIFE 20 GAUGE LANCE
MVR 21A	MVR KNIFE 21 GAUGE ANGLED
MVR 21	MVR KNIFE 21 GAUGE LANCE
MVR 23A	MVR KNIFE 23 GAUGE ANGLED
MVR 23	MVR KNIFE 23 GAUGE LANCE
MVR 25A	MVR KNIFE 25 GAUGE ANGLED
MVR 25	MVR KNIFE 25 GAUGE LANCE

## Handling Instructions



**Taking out the knife:** Take special care to prevent the blade from coming into contact with any part of the packaging. Gently pick up the device and lift the knife from the locking position.

**Replacing the knife:** Hold the knife with the blade edge facing downward. Align it properly with the locking position in the tray. Press down gently until the handle clicks into place, ensuring the knife is securely locked.

## Order Samples

Visit our website to request a sample of our surgical instruments and experience firsthand the commitment to quality that defines Makvin Surgical.



## Contact Us



## Makvin Surgical LLC

Phone : +1 423 750 4247

E-mail : [sales@makvinsurgical.com](mailto:sales@makvinsurgical.com)

Website : [www.makvinsurgical.com](http://www.makvinsurgical.com)

# MAKVIN Surgical